

# Quick Guide

## MYgen Gel & PCR Purification System - Agarose Gel Extraction

Products	Cat No.	Size
MYgen Gel & PCR Purification System	MYG104-100	100 preps
	MYG104-200	200 preps

### Kit Contents:

UB  
WB Bottle  
Help B Buffer  
EB DNA Hydration Solution  
Spin column / Collection tube

### Preparation:

- Prepare **80% Ethanol** (not provided) in the **WB** bottle.
- Prepare **Isopropanol** (not provided)


*Note: Make fresh before use.*

### Protocol:

1. Weigh a new empty 1.5 mL tube. Then excise or cut out the selected DNA fragment of your choice from the agarose gel
  - Transfer the excised DNA gel fragment to the weighed 1.5 mL tube and weigh the gel fragment.
  - Add 3 volumes of **UB** buffer to 1 volume of gel  
*\*Eg. Add 300 uL of UB buffer to 100 mg of gel fragment*
  - Incubate the mixture for 10 mins at 65°C
  - Invert the mixture during incubation to fully dissolve gel
  - Add 1 volume of **Isopropanol** and mix well  
*\*Eg. Add 100 uL of Isopropanol to 100 mg of gel fragment*

*Note: For >2% gel, kindly use 6 volumes of UB buffer*
2. Place the spin column in the 2 mL collection tube provided
  - Add 200 uL of **Help B** buffer to spin column
  - Centrifugation at 10,000 rpm for 30 seconds
  - Discard flow-through, place back the spin column into the 2 mL collection tube
3. Transfer mixture from **Step 1** into the spin column with 2 mL collection tube from **Step 2**
  - Centrifugation at 7,000 rpm for 1 minute
  - Discard flow-through, place back the spin column into the 2 mL collection tube

*\* To obtain high yield DNA recovery, add flow-through back in spin column and repeat step 3*
4. Add 750 uL **WB** (80% Ethanol) to the spin column
  - Centrifuge at 13,000 rpm for 30 sec
  - Discard flow-through, place back the spin column into 2 mL collection tube
  - Repeat Step 4 washing again
5. Centrifuge at 13,000 rpm for 3 mins to remove leftover residues
6. Place the spin column in a new 1.5 mL micro tube (not provided)
  - Discard the used 2 mL collection tube
7. Add 30 to 50 uL buffer **EB** DNA Hydration Solution
  - Incubate the mixture for 1 min at RT
  - Centrifuge at 13,000 rpm for 2 minute
  - Discard the spin column
  - Store at 4°C or -20°C



# Quick Guide

## MYgen Gel & PCR Purification System - PCR Clean Up

Products	Cat No.	Size
MYgen Gel & PCR Purification System	MYG104-100	100 preps
	MYG104-200	200 preps

### Kit Contents:

UB  
WB Bottle  
Help B Buffer  
EB DNA Hydration Solution  
Spin column / Collection tube

### Preparation:

- Prepare **80% Ethanol** (not provided) in the **WB** bottle.
- Prepare **Isopropanol** (not provided)

*Note: Make fresh before use.*

Sample Preparation: For high recovery yield of PCR product

- Use 1 volume of PCR sample (50 uL)
- Add 3 volumes of **UB** buffer (150 uL)
- Add 2 volumes of **Isopropanol** (100 uL)

### Protocol:

1. Mix PCR samples as mentioned in Sample Preparation step above in a new 1.5 mL tube (not provided)
2. Place the spin column in the 2 mL collection tube provided
  - Add 200 uL of **Help B** buffer to spin column
  - Centrifugation at 10,000 rpm for 30 seconds
  - Discard flow-through, place back the spin column into the 2 mL collection tube
3. Transfer mixture from **Step 1** into the spin column with 2 mL collection tube from **Step 2**
  - Centrifugation at 7,000 rpm for 1 minute
  - Discard flow-through, place back the spin column into the 2 mL collection tube

*\* To obtain high yield DNA recovery, add flow-through back in spin column and repeat step 3*
4. Add 750 uL **WB** (80% Ethanol) to the spin column
  - Centrifuge at 13,000 rpm for 30 seconds
  - Discard flow-through, place back the spin column into 2 mL collection tube
  - Repeat Step 4 washing again
5. Centrifuge at 13,000 rpm for 3 minutes to remove leftover residues
6. Place the spin column in a new 1.5 mL micro tube (not provided)
  - Discard the used 2 mL collection tube
7. Add 30 to 50 uL buffer **EB** DNA Hydration Solution
  - Incubate the mixture for 1 minute at room temperature
  - Centrifuge at 13,000 rpm for 2 minute
  - Discard the spin column
  - Store at 4°C or -20°C

1. Use 1 volume of PCR sample (50 uL)  
Add 3 volumes **UB** (150 uL)  
Add 2 volumes **Isopropanol** (100 uL)
2. Place spin column in 2 mL collection tube  
Add 200 uL **Help B** buffer to spin column  
Centrifuge (10,000 rpm, 30 sec)  
Discard flow through
3. Transfer **Step 1** mixture into spin column from **Step 2**  
Centrifuge (7,000 rpm, 1 min)  
Discard flow through
4. Add 750 uL **WB** into spin column  
Centrifuge (13,000 rpm, 30 sec)  
Discard flow-through  
Repeat WB washing step again
5. Centrifuge (13,000 rpm, 3 mins)
6. Place spin column into new 1.5 mL tube  
Discard used 2 mL collection tube
7. Add 30-50 uL Buffer **EB**  
Incubate mixture (RT, 1 min)  
Centrifuge (13,000 rpm, 2 min)

DNA Elution  
Store at 4°C or -20°C



# Quick Guide

## MYgen Gel & PCR Purification System - Primer Dimer Removal

Products	Cat No.	Size
MYgen Gel & PCR Purification System	MYG104-100	100 preps
	MYG104-200	200 preps

### Kit Contents:

UB  
WB Bottle  
Help B Buffer  
EB DNA Hydration Solution  
Spin column / Collection tube

### Preparation:

- Prepare **80% Ethanol** (not provided) in the **WB** bottle.

*Note: Make fresh before use.*

Sample Preparation: For dimer removal from PCR product

- Use 1 volume of PCR sample (50 uL)
- Add 5 volumes of **UB** buffer (250 uL)

### Protocol:

1. Mix PCR samples as mentioned in Sample Preparation step above in a new 1.5 mL tube (not provided)
2. Place the spin column in the 2 mL collection tube provided
3. Transfer mixture from **Step 1** into the spin column with 2 mL collection tube from **Step 2**
  - Centrifugation at 7,000 rpm for 1 minute
  - Discard flow-through, place back the spin column into the 2 mL collection tube

*\* To obtain high yield DNA recovery, add flow-through back in spin column and repeat step 3*
4. Add 750 uL **WB** (80% Ethanol) to the spin column
  - Centrifuge at 13,000 rpm for 30 seconds
  - Discard flow-through, place back the spin column into 2 mL collection tube
  - Repeat Step 4 washing again
5. Centrifuge at 13,000 rpm for 3 minutes to remove leftover residues
6. Place the spin column in a new 1.5 mL micro tube (not provided)
  - Discard the used 2 mL collection tube
7. Add 30 to 50 uL buffer **EB** DNA Hydration Solution
  - Incubate the mixture for 1 minute at room temperature
  - Centrifuge at 13,000 rpm for 2 minute
  - Discard the spin column
  - Store at 4°C or -20°C

