

Quick Guide

MYgen Plasmid Mini Prep Kit

Products	Cat No.	Size
MYgen Plasmid Mini Prep Kit	MYP101-100	100 preps
	MYP101-200	200 preps

Kit Contents:

B1
B2
RNase A
Blue Indicator
WB Bottle
EB DNA Hydration Solution
Spin column / Collection tube

Preparation:

- Prepare 80% Ethanol (not provided) in the WB bottle.
Note: Make it fresh directly before use.
 - Add Deionized Water (D. W) to RNase A tube according to label. Transfer mixed RNase A solution to B2 and mix well before use. Store at 4°C after mixed.
 - Add Ethanol to Blue Indicator tube according to label. Transfer mixed solution to B1 and mix well before use. (Optional)
- * Note: Blue indicator is a pH indicator which appears blue color. This indicator does not affect the extraction efficiency or purity.

Protocol:

1. Use 1 to 3 mL microbial culture
 - Centrifuge at 10,000 rpm for 1 min
 - Leave about 100ul of media + pellet and discard the rest of the supernatant
 - Add 350 ul of B1 (Optional: Blue Indicator) to cell pellet.
 - Mix well by vortexing for 5 seconds.
2. Add 350 ul B2 to the mixture from Step 1
 - Mix gently by inverting 4 to 6 times (Caution: Do not vortex)
 - Centrifuge at 13,000 rpm for 3 to 5 mins or leave on ice for 5 mins
3. Transfer supernatant from Step 2 into the spin column with 2 mL collection tube that is provided
 - Centrifugation at 7,000 rpm for 30 seconds
 - Discard flow-through, place back the spin column into the 2 mL collection tube
4. Add 750 ul WB (80% Ethanol added) to the spin column
 - Centrifuge at 13,000 rpm for 30 sec
 - Discard flow-through, place back the spin column into 2 mL collection tube
 - Repeat Step 4 washing again
5. Centrifuge again at 13,000 rpm for 3 mins to remove leftover residues
6. Place the spin column in a new 1.5 mL micro tube (not provided)
 - Discard the used 2 mL collection tube
7. Add 30 to 50 ul Buffer EB DNA Hydration Solution
 - Incubate the mixture for 1 min at RT
 - Centrifuge at 13,000 rpm for 1 min
 - Discard the spin column
 - Store at 4°C or -20°C

