

# Quick Guide

## MYgen Plasmid Mini Prep Kit

Products	Cat No.	Size
MYgen Plasmid Mini Prep Kit	MYP101-100	100 preps
	MYP101-200	200 preps

### Kit Contents:

B1	Help B Buffer
B2	WB Bottle
RNase A	EB DNA Hydration Solution
Blue Indicator	Spin column / Collection tube

### Preparation:

- Prepare **80% Ethanol** (not provided) in the **WB bottle**.  
Note: Make it fresh directly before use.
- Add Deionized Water (D. W) to **RNase A** tube according to label. Transfer mixed RNase A solution to **B2** and mix well before use. **Store at 4°C after mixed.**
- Add **Ethanol** to **Blue Indicator** tube according to label. Transfer mixed solution to **B1** and mix well before use. (Optional)  
*\* Note: Blue indicator is a pH indicator which appears blue color. This indicator does not affect the extraction efficiency or purity.*

### Protocol:

1. Use 1 to 3 mL microbial culture
  - Centrifuge at 10,000 rpm for 1 min
  - Leave about 100 ul to 200 ul of media + pellet and discard the rest of the supernatant.
  - Resuspend pellet by vortexing for 10-30 seconds.
2. Add 350 ul of **B1** (Optional: Blue Indicator) to cell pellet.
  - Mix well by inverting 5 to 10 times.
  - Add 350 ul **B2** to the mixture from Step 1
  - Mix gently by inverting 5 to 10 times (**Caution: Do not vortex**)
  - Centrifuge at 1,000 rpm for 5 mins.
3. Separately place the **spin column** into **2 ml collection tube** provided.
  - Add 200 ul **Help B** buffer to spin column.
  - Centrifuge at 10,000 rpm for 30 seconds.
  - Discard flow through.
4. Transfer **supernatant from Step 2** into the spin column with 2 mL collection tube from Step 3.
  - Centrifugation at 7,000 rpm for 1 min
  - \* Note: Repeat step 4 with flow-through solution to obtain higher yield plasmid*
  - Discard flow-through, place back the spin column into the 2 mL collection tube
5. Add 750 ul **WB** (80% Ethanol added) to the spin column
  - Centrifuge at 13,000 rpm for 30 sec
  - Discard flow-through
  - Repeat Step 5 washing again
6. Centrifuge again at 13,000 rpm for 3 mins to remove leftover residues
7. Place the spin column in a new 1.5 mL micro tube (not provided)
  - Discard the used 2 mL collection tube
8. Add 30 to 50 ul Buffer EB DNA Hydration Solution
  - Incubate the mixture for 1 min at RT
  - Centrifuge at 13,000 rpm for 1 min
  - Discard the spin column
  - Store at 4°C or -20°C

