


Quick Guide

MYgen Total RNA Prep Kit - Tissue Culture Cell Protocol

Products	Cat No.	Size
MYgen Total RNA Prep Kit (For Bacterium, Plant, Animal Tissue, Tissue cultured Cell)	MYR101-050	50 preps
	MYR101-100	100 preps

Kit Contents:

RB, Proteinase K
RW, RNase Free Water
HiSpin Capsule



Preparation:

- 2-Mercaptoethanol(2-ME) should be added to RB solution and the mixed RB/2-ME Solution should be stored at 4°C. (Ratio of addition: 10 μ L 2-ME : 1mL RB Solution)
- Add 100% ethanol to RW according to label before use.
- Add Deionized Water to lyophilised Proteinase K according to label and mix well. Store at -20°C.

Protocol:

1. Vortex mixture of harvest cells ($< 1 \times 10^7$) + 350 μ L RB(+2-ME) + 5 μ L Proteinase K and incubation at 56°C for 10 mins.
2. Perform centrifugation at 14,000 rpm for 3 mins, 4°C and transfer the supernatant into a new 1.5mL micro tube.
3. Add 250 μ L Ethanol (100%) into lysate and vortex for 30 sec – 1 mins
4. Place the spin column into a 2 mL collection tube then
 - Add Lysate/Ethanol mix from Step 3 into spin column
 - Perform centrifugation at 14,000 rpm for 30 sec, 4°C
 - Discard the flow-through and place the spin column into a 2 mL collection tube again
5. Add 500 μ L RW (RNA washing solution) into spin column
 - Perform centrifugation at 14,000 rpm for 30 sec, 4°C
 - Discard the flow-through and place the spin column into a 2 mL collection tube again
 - Repeat above steps for washing
6. Perform centrifugation at 14,000 rpm for 3 mins, 4°C to remove residual wash buffer
 - Place the spin column into a new 1.5 mL micro tube
7. Add 50 - 100 μ L RNase free water
 - Incubate the mixture for 1 min at RT
 - Perform centrifugation at 14,000 rpm for 1 min
 - Remove the Spin column
 - Store at -70°C




Quick Guide

MYgen Total RNA Prep Kit - Plant Tissue Protocol

Products	Cat No.	Size
MYgen Total RNA Prep Kit (For Bacterium, Plant, Animal Tissue, Tissue cultured Cell)	MYR101-050	50 preps
	MYR101-100	100 preps

Kit Contents:

RB, Proteinase K
RW, RNase Free Water
HiSpin Capsule



Preparation:

- 2-Mercaptoethanol(2-ME) should be added to RB solution and the mixed RB/2-ME Solution should be stored at 4°C. (Ratio of addition: 10 μ L 2-ME : 1mL RB Solution)
- Add 100% ethanol to RW before use.
- Sample
(It is recommended to start with no more than 100 mg plant tissue. Grind the fresh-freeze tissues with liquid nitrogen thoroughly. Keep the samples cool and perform the following steps promptly to minimize RNA degradation)

Protocol:

1. Vortex mixture sample (< 100 mg) + RB(+2-ME) 350 μ L for 30 sec - 1 min and incubation at 60°C for 3 mins
(Note: Perform incubation at 60°C for 5 - 10 mins in case of high starch or polyphenol / polysaccharide sample)
2. Perform centrifugation at 14,000 rpm for 3 mins, 4°C and transfer the supernatant into a new 1.5mL micro tube.
3. Add 250 μ L Ethanol (100%) into lysate and vortex for 30 sec - 1 mins
4. Place the spin column into a 2 mL collection tube then
- Add Lysate/Ethanol mix from Step 3 into spin column
- Perform centrifugation at 14,000 rpm for 30 sec, 4°C
- Discard the flow-through and place the spin column into a 2 mL collection tube again
5. Add 500 μ L RW (RNA washing solution) into spin column
- Perform centrifugation at 14,000 rpm for 30 sec, 4°C
- Discard the flow-through and place the spin column into a 2 mL collection tube again
- Repeat above steps for washing
6. Perform centrifugation at 14,000 rpm for 3 mins, 4°C to remove residual wash buffer
- Place the spin column into a new 1.5 mL micro tube
7. Add 50 - 100 μ L RNase free water
- Incubate the mixture for 1 min at RT
- Perform centrifugation at 14,000 rpm for 1 min
- Remove the Spin column
- Store at -70°C



Quick Guide

MYgen Total RNA Prep Kit - Bacterium Protocol

Products	Cat No.	Size
MYgen Total RNA Prep Kit (For Bacterium, Plant, Animal Tissue, Tissue cultured Cell)	MYR101-050	50 preps
	MYR101-100	100 preps

Kit Contents:

RB, Proteinase K
RW, RNase Free Water
HiSpin Capsule

1. Harvest Bacteria ($< 1 \times 10^9$)
Add 100 uL TE
Add 2 uL of Lysozyme
Vortex (30 sec - 1 min)
Incubation (RT, 3 - 10 mins)
2. Add 350 uL RB1 (+2-ME)
Vortex (30 sec - 10 mins)
Cfg (14,000 rpm, 3 mins)
Transfer supernatant to new 1.5mL micro tube
3. Add 250 uL Ethanol (100%)
Vortex (30 sec - 1 min)
4. Add lysate/ethanol mix to spin column
Cfg (14,000 rpm, 30 sec, 4°C)
Discard flow-through
5. Add 500 uL RW into spin column
Cfg (14,000 rpm, 30 sec, 4°C)
Discard flow-through
Repeat RW washing step again
6. Cfg (14,000 rpm, 3 mins, 4°C)
- remove final wash buffer residual
Replace a new 1.5mL micro tube
7. Add 50-100 uL RNase free water
Incubate mixture (RT, 1 min)
Cfg (14,000 rpm, 1min)

RNA Elution
Store at -70°C

Preparation:

- 2-Mercaptoethanol(2-ME) should be added to RB solution and the mixed RB/2-ME Solution should be stored at 4°C. (Ratio of addition: 10 uL 2-ME : 1mL RB Solution)
- Add 100% ethanol to RW according to label before use.
- Add Deionized Water to lyophilised Lysozyme according to label and mix well. Store at -20°C.
- Sample
 - i. Gram-negative bacteria: ~ 500 ug
 - ii. Gram-positive bacteria: ~ 3 mg / mL

Protocol:

1. Vortex mixture of harvest bacteria ($< 1 \times 10^9$) + TE buffer 100 uL + Lysozyme 2 uL and incubate:
 - i. Gram-negative bacteria: 3 - 5 mins, RT
 - ii. Gram-positive bacteria: 5 - 10 mins, RT
2. Add in 350 uL RB(+2-ME) and vortex for 30 sec - 1 min
- Perform centrifugation at 14,000 rpm for 3 mins, 4°C and transfer the supernatant into a new 1.5mL micro tube.
3. Add 250 uL Ethanol (100%) into lysate and vortex for 30 sec - 1 mins
4. Place the spin column into a 2 mL collection tube then
- Add Lysate/Ethanol mix from Step 3 into spin column
- Perform centrifugation at 14,000 rpm for 30 sec, 4°C
- Discard the flow-through and place the spin column into a 2 mL collection tube again
5. Add 500 uL RW (RNA washing solution) into spin column
- Perform centrifugation at 14,000 rpm for 30 sec, 4°C
- Discard the flow-through and place the spin column into a 2 mL collection tube again
- Repeat above steps for washing
6. Perform centrifugation at 14,000 rpm for 3 mins, 4°C to remove residual wash buffer
- Place the spin column into a new 1.5 mL micro tube
7. Add 50 - 100 uL RNase free water
- Incubate the mixture for 1 min at RT
- Perform centrifugation at 14,000 rpm for 1 min
- Remove the Spin column
- Store at -70°C




Quick Guide

MYgen Total RNA Prep Kit - Animal Tissue Protocol

Products	Cat No.	Size
MYgen Total RNA Prep Kit (For Bacterium, Plant, Animal Tissue, Tissue cultured Cell)	MYR101-050	50 preps
	MYR101-100	100 preps

Kit Contents:

RB, Proteinase K
RW, RNase Free Water
HiSpin Capsule



Preparation:

- 2-Mercaptoethanol(2-ME) should be added to RB solution and the mixed RB/2-ME Solution should be stored at 4°C. (Ratio of addition: 10 uL 2-ME : 1mL RB Solution)
- Add 100% ethanol to RW according to label before use.
- Add Deionized Water to lyophilised Proteinase K according to label and mix well. Store at -20°C.
- Sample
(It is recommended to start with no more than 30 mg animal tissue. Grind the fresh-freeze tissues with liquid nitrogen thoroughly. Keep the samples cool and perform the following steps promptly to minimize RNA degradation)

Protocol:

1. Vortex mixture sample (< 30 mg) + RB(+2-ME) 350 uL + Pro-K (10 mg/ml) 5uL for 30 sec - 1 min and incubation at 56°C for 10 mins
2. Perform centrifugation at 14,000 rpm for 3 mins, 4°C and transfer the supernatant into a new 1.5mL micro tube.
3. Add 250 uL Ethanol (100%) into lysate and vortex for 30 sec – 1 mins
4. Place the spin column into a 2 mL collection tube then
 - Add Lysate/Ethanol mix from Step 3 into spin column
 - Perform centrifugation at 14,000 rpm for 30 sec, 4°C
 - Discard the flow-through and place the spin column into a 2 mL collection tube again
5. Add 500 uL RW (RNA washing solution) into spin column
 - Perform centrifugation at 14,000 rpm for 30 sec, 4°C
 - Discard the flow-through and place the spin column into a 2 mL collection tube again
 - Repeat above steps for washing
6. Perform centrifugation at 14,000 rpm for 3 mins, 4°C to remove residual wash buffer
 - Place the spin column into a new 1.5 mL micro tube
7. Add 50 - 100 uL RNase free water
 - Incubate the mixture for 1 min at RT
 - Perform centrifugation at 14,000 rpm for 1 min
 - Remove the Spin column
 - Store at -70°C

